

MODEL
YE246ES

ULTRA – ECONOMIC MODEL



INTELLIGENT CONTROL FOR 360°



1. Outdoor and indoor using.
2. Support heavy loading capacity 135KGS.
3. Intelligent control for 360° rotation.
4. Outstanding matching color makes it attractive.
5. Joystick on left or right side.



Product parameters



Model No.	Recline
YE246ES	NO
Seat width	Seat depth
45 cm	48 cm
Front wheel	Rear wheel
8" PU Solid wheel	12" PU Solid wheel
Battery	Motor
Lithium Battery 11Ah	Brushed 250W each
Charger	Controller
24V 2A	Brushed controller
Brakes	Drive distance
Electromagnetic braking	12-15 km
Max speed	Net weight
6 km/h	25 kg
Weight capacity	Unfolded Size
135 kg	105*60*94cm

EASY CLIMBING

