



# ECONOMIC FOLDABLE WHEELCHAIR



- ELECTROMAGNETIC BRAKE
- MANUAL BRAKE
- STRONG STRRL FRAM
- 360° JOYSTICK CONTROLLER

[www.yattli.com](http://www.yattli.com)

# UPDATED DEVICE WITHOUT RISE IN PRICE

This wheelchair is made of strong steel material. The classic black color with silver rim, makes the wheelchair classic.

This chair has flip up armrest, foldable legrest and removable cushions, which makes the maintenance simple for daily life.

It comes with 24V 24Ah battery which lasting for 18KM, and the special price makes it definitely selling good in Asia countries and Mid-east countries.



Foldable backrest



Flip up armrest



High endurance batteries



Foldable legrest



Anti-tip



High strength fork

**YATTLI INDUSTRY CO., LTD.**



Tel: +86 755 817 071 99

Web: [www.yattli.com](http://www.yattli.com)

Email: [info@yattli.com](mailto:info@yattli.com)

Factory: No.2 Lingchang Road, PingShan Village, TangXia Town, DongGuang City, GuangDong, China

Office: RM.1403, BLK.B of Duocai Office Complex, GuanlanHigh-tech Zone, Longhua DIST., Shenzhen, China

# INTERLACED MOTORS DESIGN



## Features

- Strong steel frame
- Programmable controller, 360° control
- Flip-up armrest and foldable legrest
- With 11AH lead acid battery
- Electromagnetic braking system
- With 12.5" rear wheel, strong with aluminum rim

## Parameter

**Front wheel:** 8", PU Solid Tyre

**Rear Wheel:** 12.5", Aluminum Rim

**Net Weight:** 42 kg

**Brake:** Electromagnetic Brake

**Batteries:** Lead Acid 24V 11Ah

**Motors (pair):** Brushed Motor 250W

**Charger:** 24V 2A

**Controller:** Programmable Controller

**Maximum Range:** 15-18 km

**Ground Clearance:** 8 cm

**Max Speed:** 6 km/h

**Charging Hour:** 6-8 Hours

**Max Slope:** 12 degrees

**Weight Capacity:** 125 KG

**Dimension (Folded):** 79\*72\*37cm



## Footnotes

1) Range and speed vary with userweight-terrain type, battery charge battery condition and tyre pressure.

2) Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.

3) Battery weight may vary based on manufacturer.

**YATTLL INDUSTRY CO., LTD.**

Tel: +86 755 817 071 99

Web: [www.yattll.com](http://www.yattll.com)

Email: [info@yattll.com](mailto:info@yattll.com)

Factory: No.2 Lingchang Road, PingShan Village, TangXia Town, DongGuang City, GuangDong, China

Office: RM.1403, BLK.B of Duocai Office Complex, GuanlanHigh-tech Zone, Longhua DIST., Shenzhen, China

