


Reinforced Frame

Loading capacity 150KG



Main Features:

1. Redesigned and reinforced frame which can load 150KG.
2. Powerful & brushless motors 250W*2.
3. With top advanced seat back cushion.
4. It has a quick-reased battery.
5. 8" rubber wheel, 12.5" rear wheel.
6. Joystick can be mounted left and right side.




Reinforced Frame

Stronger and safer



360°

Intelligent controller



Suspension System

Double shock absorber spring



Perforated Wheels

With absorbing ability



Product parameters



Model No. YE100	Recline NO
Seat width 49 / 52 / 56 cm	Seat depth 50 cm
Front wheel 8" Suspension wheel	Rear wheel 12.5" Solid, Suspension wheel
Battery 11Ah/15Ah	Motor Brushless 250W
Charger 24V 2A/3A/5A	Controller Brushless controller
Brakes Electromagnetic braking	Drive distance 15-18 km
Max speed 6 km/h	Net weight 24.5 kg
Weight capacity 150 kg	Expanded Size 108*60*99 cm

Email: info@yattll.com

Factory: YALU TECHNOLOGY PARK, NO. 88 Pingshan Industrial Rd, Tangxia, Dongguan City, Guangdong, China.

Office: Rm 907, Yichuang International Centre, No. 8 Gaoerfu blvd, Longhua Dist, Shenzhen, Guangdong, China



YATLL INDUSTRY CO.,LTD

755 - 817 - 07199 www.yattll.com

