

CHOOSE THE COLOUR THAT SUITS YOU



REBOOT YOUR LIFE

Model No.	YE100R
Seat width	49 / 52 / 56 cm
Seat depth	44 cm
Front wheel	8" Suspension wheel
Rear wheel	12.5" Suspension wheel
Battery	24V, 10.4Ah
Motor	Brushless motors 180W*2
Charger	24V 2A/3A/5A
Controller	Brushless controller
Brakes	Electromagnetic braking
Drive distance	15-18 km
Recline	Yes
Weight capacity	150 kg
Max speed	8 km/h
Net weight	29.6kg
Unfolded Size	108*64*110 cm

Main Features:

1. It has electric-recline function
2. Different angel-anjustable headrest
3. With medical cylinder in the back supporting the recline
4. Upgraded brushless controller, 360° rotation
5. Featuring unique wheels with peforations that has strong absorbing ability
6. Supporting loading capacity 150KG



Electric recline by remote control



Reinforced Frame Loading Capacity 150KG



Quick Folding



Economical
Model **YE100R**

Intelligent remote control

