

Electrical Recline

YE200R



Product Details



Brushless controller with recline function



Extended Footpad



Angle-Adjustable Headrest



Unique Dual-battery Design





Intelligent Control

Intelligent controller with remote control, easy to use



Brushless Motors

Powerful brushless motors 250W*2, easy travelling



Suspension System

Dual suspension system can well prevent turbulences



12 Ah / 18Ah optional



Product Parameters



Model No.

YE200R

Recline

Yes

Seat width

45 / 52 / 56 cm

Seat depth

43 cm

Front wheel

8" Solid wheel

Rear wheel

12.5" Solid, Aluminum rim

Battery

6Ah / 9Ah *2pcs

Motor

Brushless motors 250W*2

Charger

24V 2A/3Ah/5Ah

Controller

Brushless controller

Brakes

Electromagnetic braking

Drive distance

15-18 km

Max speed

6 km/h

Net weight

31.6 kg

Weight capacity

150 kg

Expanded Size

110*60*98cm