

Aluminium Frame **YEAZI**



STRONG ALUMINIUM FRAME

Powerful and Lightweight.



**One-step
Folding**



**Net Weight
- 16kg**



**Durable PU
solid wheels**



**Brushless
Controller**



**Lithium
Battery
10.4Ah**



**Brushless
Motors
180W*2**

Easy-to-Use Intelligent Controller

Battery Indicator

**Speed Down / “-”
Speed Up / “+”**

Horn

Power On/Off

Speed Level

Light



Product Parameters

Model No. YEAZ1	Recline NO
Seat width 45 cm	Seat depth 40 cm
Front wheel 7 inches	Rear wheel 8 inches
Battery Lithium 10.4Ah	Motor Brushless 180W *2
Charger 24V / 2A	Controller Brushless controller
Brakes Brushless controller	Drive distance 16-20 km
Max speed 6 km/h	Net weight 16 kg
Weight capacity 125 kg	Expanded Size 975*605*945cm