


NEW GENERATION OF YE246



Main Features:

- Fashionable & Comfortable seat back cushion
- Widened & upgraded armrest
- Features unique wheels with hundreds of perforations
- Brushless motors 250W*2
- Loading capacity 135kg
- Max driving range 15-18km



**Widened &
Comfortable armrest**



Brushless controller



**Perforated wheels
with absorbing ability**



**Brushless motors
250W*2**

Fashionable & Comfortable model

2023



Product parameters

Model No.	Recline
YE246 II	NO
Seat width	Seat depth
45.5 cm	42 cm
Front wheel	Rear wheel
8" solid wheel	12.5" solid, Suspension wheel
Battery	Motor
11Ah/15Ah/20Ah	Brushless 250W
Charger	Controller
24V 2A/3A/5A	Brushless controller
Brakes	Drive distance
Electromagnetic braking	15-18 km
Max speed	Net weight
6 km/h	25.5kg
Weight capacity	Expanded Size
135 kg	61*78*100 cm

Email: info@yattll.com

Factory: YALU TECHNOLOGY PARK, NO. 88 Pingshan Industrial Rd, Tangxia, Dongguan City, Guangdong, China.

Office: Rm 907, Yichuang International Centre, No. 8 Gaoerfu Blvd, Longhua Dist, Shenzhen, Guangdong, China



YATLL INDUSTRY CO.,LTD

755 - 817 - 07199 www.yattll.com

